

We invented the dermatoscope 30 years ago. Now we've done it again.

The new Dermatoscopes: DELTA 30 and DELTAone.



30
YEARS
1st DERMATOSCOPE

M0223 02

Distributed by
HENLEYS MEDICAL SUPPLIES LIMITED
Brownfields, Welwyn Garden City, Herts AL7 1AN
Tel: 01707 333164 Fax: 01707 334795
www.henleysmed.com



THE NEW HEINE DELTA 30 DERMATOSCOPE.

**If there's something to discover,
you'll be sure to see it.**

Not only would our new high-end dermatoscope win any beauty pageant; it stands out due to the best imaging we've ever developed. Combined with a field of view measuring a veritable 30 millimetres and unprecedented colour rendering, in LED^{HQ} of course.

Polarisation provides glare-free and non-reflective working conditions. Above all else, it is comfortable to operate with its angled, ergonomic design and intuitive control. Feel free to take a look – you're in for a real discovery.



THE NEW HEINE DELTA*one* DERMATOSCOPE.

Mobility and quality have now become *one*.

Small on the outside, big on the inside – this mobile masterpiece is taking medical professionals by storm. The crystal-clear image was accomplished using the new achromatic HEINE optical system. What's more, it produces extraordinary colour rendering, due to the quality of our LED^{HQ}. Examinations are now particularly

comfortable due to the polarisation allowing doctors to work glare-free and without disrupting reflections. These three features – clarity, LED^{HQ} and non-reflective working – are what turn the DELTA*one* into a vital tool to carry in your pocket wherever you go.

Crystal-clear image with the new achromatic HEINE optical system

Fits into any pocket

Glare-free and non-reflective working conditions thanks to polarisation

LED^{HQ}
LED NOW IN HEINE QUALITY.

Excellent colour rendering due to LED^{HQ}

Digital documentation with the exclusive HEINE DERM app



THE NEW DERMATOSCOPES: DELTA 30 AND DELTAone.

Seeing is believing.

Colours just as God created them. The DELTA 30 provides the best image quality of any dermatoscope we've ever built. Due to its small dimensions – and high diagnostic precision – the DELTAone is always a tool that users will gladly carry in their pocket. Their full


potential is revealed only by looking through them: the new achromatic HEINE optical systems provide a level of colour rendering that has never been seen before thanks to the specially developed LED^{HQ}.




Conventional LED



LED^{HQ}
LED NOW IN HEINE QUALITY.



Conventional presentation



Presentation with one of the new DELTA Dermatoscopes

When it comes to the minutest of detail.

The new achromatic HEINE optical system, installed in the two new dermatoscopes of the DELTA series, provides optimum illumination right to the outer edges. The light is distributed homogeneously and evenly across the entire

field of view. The result is no disturbing reflections, shadows or bright or dark spots. With unseen clarity and an outstanding level of detail – regardless of the working distance. This feature is available in the new DELTA 30 with an impressive field of view size of real 30mm in diameter.

THE NEW DERMATOSCOPES: DELTA 30 AND DELTAone.

Simply digital. Digitally simple.

Both of the new dermatoscopes can also be combined with an iPhone* from version 6 upwards using the (also new) adapter. Our exclusive HEINE DERM App lets you create patient files, manage and archive pictures,

and save the serial documentation of progression over time. The HEINE DERM App is easy to install and the software is simple and intuitive to use.



*The Apple iPhone is not included in the scope of delivery. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. **Not available in all countries.

ORDERING

Let's get numerical.

DELTA 30 with contact plate with scale, USB cord with medical approved plug-in power supply, case	K-230.28.305
DELTAone Standard Edition (black-silver) with contact plate with scale, USB cord with medical approved plug-in power supply, case	K-210.28.305
DELTAone White Edition (white-silver) with contact plate with scale, USB cord with medical approved plug-in power supply, case	K-211.28.305
Contact plate with scale for DELTA 30	K-000.34.220
Small contact plate for DELTA 30. For the examination of difficult-to-access pigmented lesions.	K-000.34.206
Contact plate with scale for DELTAone	K-000.34.222
Small contact plate for DELTAone. For the examination of difficult-to-access pigmented lesions.	K-000.24.207
Table charger for DELTA 30	X-002.99.212
Mounting case smartphone for Apple iPhone 6/6s* suitable for DELTA 30/DELTAone	K-000.34.260
Mounting case smartphone for Apple iPhone 7/8* suitable for DELTA 30/DELTAone	K-000.34.261
Mounting case smartphone for Apple iPhone X/XS* suitable for DELTA 30/DELTAone	K-000.34.262
Mounting case smartphone for Apple iPhone XR* suitable for DELTA 30/DELTAone	K-000.34.263
Universal Smartphone Connector (Stick-On) . Suitable for many common smartphone mounting cases. Connector suitable for mobile phone cases with a smooth surface where the camera mustn't protrude from the case.	K-000.34.270
DELTA 30 Universal Smartphone Adaptor (Clamp-on) . Suitable for many common smartphones where a removable clamp is preferred	K-000.34.272
DELTAone Universal Smartphone Adaptor (Clamp-on) . Suitable for many common smartphones where a removable clamp is preferred	K-000.34.271

*The Apple iPhone is not included in the scope of delivery. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. **Not available in all countries.

ORDER ON-LINE

www.henleysmedicalsupplies.co.uk/dermatoscopes



M0223 02

Distributed by
HENLEYS MEDICAL SUPPLIES LIMITED
 Brownfields, Welwyn Garden City, Herts AL7 1AN
 Tel: 01707 333164 Fax: 01707 334795
www.henleysmed.com



10% off HEINE Dermatoscopes

Use code **10DELTA**

Enter the code at

www.henleysmedicalsupplies.co.uk

Offer can only be used once per customer and applies to HEINE DELTAone and HEINE DELTA30 purchased on-line only. This offer maybe subject to change or termination at any time without prior notification.