

HEINE EL 10 LED® Examination Light

A real highlight in all respects with LED^{HQ} illumination

LED^{HQ}
LED NOW IN HEINE QUALITY



The EL10 LED examination light distinguishes itself through an absolutely homogeneous field of light. The light intensity can be regulated continuously. The diameter of the spot size can be freely adjusted between 8.5 and 14 cm, with a working distance of 30 cm. Controls with ergonomic grip allow for intuitive and fast adjustment of brightness and spot size. The compact illumination head (Ø approx. 60 mm) allows an almost coaxial illumination, particularly in difficult application situations.

- ∴ **Bright and homogeneous.** Perfect illumination in every examination situation. unobstructed view of the ear canal and tympanum.
- ∴ HEINE has set a new standard which maintains that only the best is good enough, from the selection of materials to processing, from light intensity to dimmability, and from thermal management to a colour rendering index as high as possible.
That's what we call LED in HEINE quality – or LED^{HQ}.
- ∴ **3 different mounting options.** Wall mounting, clamp for table-top mounting or on a wheeled stand.
- ∴ **Solid construction.** Reliable operation and illumination power for many years.
- ∴ **Hygiene.** Easy and efficient cleaning and disinfection.
- ∴ **Coiled goose neck.** Easy and intuitive adjustment.
- ∴ **An electrical control system** monitors all functions and guarantees safety.


∴ Continuous adjustment

∴ One-hand operation

∴ Compact design

EL 10 LED Examination Light	LED
with wall mount [01]	J-008.27.001
with clamp for table-top mounting [02]	J-008.27.002
with wheeled stand, metal base [03]	J-008.27.003

Clamp for table-top mounting for EL 10 LED	J-000.27.013
Wheeled stand, metal base with mount for EL 10 LED (Height 80 cm)	Y-096.50.010
Mount for wheeled stand [04]	J-000.27.014



230 V is the standard mains voltage. If a different voltage is required, e. g. 120V, please specify when ordering.

HEINE EL3 LED Examination Light

Simply bright LED illumination in LED^{HQ}

LED^{HQ}
LED NOW IN HEINE QUALITY.



- ✚ **Simply bright.** The EL3 LED provides bright, shadow-free illumination precisely at the point of interest.
- ✚ HEINE has set a new standard which maintains that only the best is good enough, from the selection of materials to processing, from light intensity and thermal management to a colour rendering index as high as possible.
That's what we call LED in HEINE quality – or LED^{HQ}.
- ✚ **Shadow-free illumination.** The compact illumination head (Ø approx. 60 mm) allows convenient, precise positioning, particularly in difficult application situations. Optimal positioning eliminates unwanted shadows in the illuminated field.
- ✚ **Maintenance-free and easy to clean.** Cool-Touch closed head design makes disinfection safe and easy.
- ✚ **Optimal positioning.** The light weight head and stable goose neck allow easy and intuitive adjustment and keep the light exactly in the required position.
- ✚ **Various mounting options.** Wall mounting, clamp for table-top mounting or on a wheeled stand.

- ✚ Shadow-free illumination
- ✚ LED in HQ
- ✚ Maintenance-free and easy to clean

EL3 LED Examination Light	LED
with wall mount [01]	J-008.27.011
with clamp for table-top mounting [02]	J-008.27.013
with wheeled stand, metal base [03]	J-008.27.014
Clamp for table-top mounting for EL3 LED	J-000.27.013
Wheeled stand, metal base with mount for EL3 LED (Height 80 cm)	Y-096.50.010
Mount for wheeled stand [04]	J-000.27.014

